

# PCR Troubleshooting Guide

The following guide can be used to troubleshoot PCR reactions. Use our [T<sub>m</sub> calculator](#) to help plan experiments and [click here for optimization tips](#).

Observation	Possible Cause	Solution
SEQUENCE ERRORS	Low fidelity polymerase	<ul style="list-style-type: none"> <li>Choose a higher fidelity polymerase such as Q5® (<a href="#">NEB #M0491</a>), Phusion® (<a href="#">NEB #M0530</a>) DNA Polymerases</li> </ul>
	Suboptimal reaction conditions	<ul style="list-style-type: none"> <li>Reduce number of cycles</li> <li>Decrease extension time</li> <li>Decrease Mg<sup>++</sup> concentration in the reaction</li> </ul>
	Unbalanced nucleotide concentrations	<ul style="list-style-type: none"> <li>Prepare fresh deoxynucleotide mixes</li> </ul>
	Template DNA has been damaged	<ul style="list-style-type: none"> <li>Start with a fresh template</li> <li>Try repairing DNA template with the PreCR® Repair Mix (<a href="#">NEB #M0309</a>)</li> <li>Limit UV exposure time when analyzing or excising PCR product from the gel</li> </ul>
	Desired sequence may be toxic to host	<ul style="list-style-type: none"> <li>Clone into a non-expression vector</li> <li>Use a low-copy number cloning vector</li> </ul>
INCORRECT PRODUCT SIZE	Incorrect annealing temperature	<ul style="list-style-type: none"> <li>Recalculate primer T<sub>m</sub> values using the <a href="#">NEB T<sub>m</sub> calculator</a></li> </ul>
	Mispriming	<ul style="list-style-type: none"> <li>Verify that primers have no additional complementary regions within the template DNA</li> </ul>
	Improper Mg <sup>++</sup> concentration	<ul style="list-style-type: none"> <li>Adjust Mg<sup>++</sup> concentration in 0.2–1 mM increments</li> </ul>
NO PRODUCT	Nuclease contamination	<ul style="list-style-type: none"> <li>Repeat reactions using fresh solutions</li> </ul>
	Incorrect annealing temperature	<ul style="list-style-type: none"> <li>Recalculate primer T<sub>m</sub> values using the <a href="#">NEB T<sub>m</sub> calculator</a></li> <li>Test an annealing temperature gradient, starting at 5°C below the lower T<sub>m</sub> of the primer pair</li> </ul>
	Poor primer design	<ul style="list-style-type: none"> <li>Check specific product literature for recommended primer design</li> </ul>

Observation	Possible Cause	Solution
		<ul style="list-style-type: none"> <li>• Verify that primers are non-complementary, both internally and to each other</li> <li>• Increase length of primer</li> </ul>
	Poor primer specificity	<ul style="list-style-type: none"> <li>• Verify that oligos are complementary to proper target sequence</li> </ul>
	Insufficient primer concentration	<ul style="list-style-type: none"> <li>• Primer concentration can range from 0.05–1 <math>\mu\text{M}</math> in the reaction. Please see specific product literature for ideal conditions</li> </ul>
	Missing reaction component	<ul style="list-style-type: none"> <li>• Repeat reaction setup</li> </ul>
	Suboptimal reaction conditions	<ul style="list-style-type: none"> <li>• Optimize <math>\text{Mg}^{++}</math> concentration by testing 0.2–1 mM increments</li> <li>• Thoroughly mix <math>\text{Mg}^{++}</math> solution and buffer prior to adding to the reaction</li> <li>• Optimize annealing temperature by testing an annealing temperature gradient, starting at 5°C below the lower <math>T_m</math> of the primer pair</li> </ul>
	Poor template quality	<ul style="list-style-type: none"> <li>• Analyze DNA via gel electrophoresis before and after incubation with <math>\text{Mg}^{++}</math></li> <li>• Check 260/280 ratio of DNA template</li> </ul>
	Presence of inhibitor in reaction	<ul style="list-style-type: none"> <li>• Further purify starting template by alcohol precipitation, drop dialysis or commercial clean up kit</li> <li>• Decrease sample volume</li> </ul>
	Insufficient number of cycles	<ul style="list-style-type: none"> <li>• Rerun the reaction with more cycles</li> </ul>
	Incorrect thermocycler programming	<ul style="list-style-type: none"> <li>• Check program, verify times and temperatures</li> </ul>
	Inconsistent block temperature	<ul style="list-style-type: none"> <li>• Test calibration of heating block</li> </ul>
	Contamination of reaction tubes or solutions	<ul style="list-style-type: none"> <li>• Autoclave empty reaction tubes prior to use to eliminate biological inhibitors</li> <li>• Prepare fresh solutions or use new reagents and new tubes</li> </ul>

Observation	Possible Cause	Solution
<b>MULTIPLE OR NON- SPECIFIC PRODUCTS</b>	Complex template	<ul style="list-style-type: none"> <li>• Use Q5 High-Fidelity (<a href="#">NEB #M0491</a>) or OneTaq® DNA Polymerases (<a href="#">NEB #M0480</a>)</li> <li>• For GC-rich templates, use Q5 High-Fidelity (<a href="#">NEB #M0491</a>) or OneTaq® DNA Polymerases. Include the appropriate GC enhancer.</li> <li>• For longer templates, we recommend <a href="#">LongAmp® Taq DNA Polymerase</a> or Q5 high-Fidelity DNA polymerase or Q5 Hot-Start High-Fidelity DNA Polymerase (<a href="#">NEB #M0493</a>)</li> </ul>
	Premature replication	<ul style="list-style-type: none"> <li>• Use a hot start polymerase, such as OneTaq Hot Start DNA Polymerase</li> <li>• Set up reactions on ice using chilled components and add samples to thermocycler preheated to the denaturation temperature</li> </ul>
	Primer annealing temperature too low Incorrect Mg <sup>++</sup> concentration	<ul style="list-style-type: none"> <li>• Increase annealing temperature</li> <li>• Adjust Mg<sup>++</sup> in 0.2–1 mM increments</li> </ul>
	Poor primer design	<ul style="list-style-type: none"> <li>• Check specific product literature for recommended primer design</li> <li>• Verify that primers are non-complementary, both internally and to each other</li> <li>• Increase length of primer</li> <li>• Avoid GC-rich 3' ends</li> </ul>
	Excess primer	<ul style="list-style-type: none"> <li>• Primer concentration can range from 0.05–1 μM in the reaction. Please see specific product literature for ideal conditions.</li> </ul>
	Contamination with exogenous DNA	<ul style="list-style-type: none"> <li>• Use positive displacement pipettes or non-aerosol tips</li> <li>• Set-up dedicated work area and pipettor for reaction setup</li> <li>• Wear gloves during reaction setup</li> </ul>
	Incorrect template concentration	<ul style="list-style-type: none"> <li>• For low complexity templates (i.e. plasmid, lambda, BAC DNA), use 1 pg–10 ng of DNA per 50 μl reaction</li> <li>• For higher complexity templates (i.e. genomic DNA), use 1 ng–1 μg of DNA per 50 μl reaction</li> </ul>

