


# Antibiotic Susceptibility Tests



# Antibiotics are:

- Natural or synthetic products that are used to kill or stop the growth of Bacteria :
  - bacteriostatic - stop growth (don't kill)
  - bactericidal – kill
- 

- **A test done to check the effectiveness of a drug against a bacterium and to select the best drug that acts against the bacterium.**

# Antibiotic Sensitivity Tests

Diffusion

Dilution

Diffusion & Dilution

Kirby-Bauer Method

Stokes Method

Tube Dilution

Agar Dilution

E-Test

Qualitative Methods

Quantitative Methods

# Tube dilution tests


A series of culture tubes are prepared, each containing a liquid medium and a different concentration of a chemotherapeutic agent. The tubes are then inoculated with the test organism and incubated for 16-20 hours at 35C. After incubation, the tubes are examined for turbidity (growth).

## **Minimum Inhibitory Concentration (MIC):**

Is the lowest concentration of chemotherapeutic agent capable of preventing growth of the test organism.


## **Minimum Bactericidal Concentration (MBC):**

Is the lowest concentration of the chemotherapeutic agent that results in no growth (turbidity) of the subcultures.







DAY 1



OVERNIGHT INCUBATION

DAY 2




CFU on subculture are determined by colony count:  
 $250 \text{ CFU} = 5 \times 10^5/\text{ml}$  organisms in original control tube

250 CFU

OVERNIGHT INCUBATION

DAY 3



CFU on subcultures made from non-turbid tubes are determined.  
 $7 \text{ CFU} = 70 \text{ CFU/ml}$  in tube containing  $8 \mu\text{g/ml}$  antibiotic.  $70 \text{ CFU/ml} > 0.1\%$  of  $5 \times 10^5 \text{ CFU/ml}$  in original tubes.

MBC =  $16 \mu\text{g/ml}$

## **2. The agar diffusion test (Kirby- Bauer test)**

This is one of the more commonly used methods of antimicrobial susceptibility testing. In this test, small filter paper disks (6 mm) impregnated with a standard amount of antibiotic are placed onto an agar plate to which bacteria have been swabbed.

The plates are incubated overnight, and the zone of inhibition of bacterial growth is used as a measure of susceptibility. Large zones of inhibition indicate that the organism is susceptible , while small or no zone of inhibition indicate resistance.

An interpretation of intermediate is given for zones which fall between the accepted cutoffs for the other interpretations. Many factors are involved in sensitivity disk testing and must be carefully controlled .these include size of inoculum, distribution of inoculum, incubation period, depth of the agar, diffusion rate of antibiotic, concentration of antibiotic in the disk, and growth rate of the bacterium .

**MATERIA  
LS**

Mueller-  
Hinton  
Agar


Antibiotic  
Disks

Turbidity  
Standard

Swabs



## **Material:**


- Mueller-Hinton agar plates
  - Antibiotics disks
  - sterile swabs
  - (4- to 6-hour) broth cultures of *Staphylococcus aureus*, *Escherichia coli*, *Pseudomonas aeruginosa*, and *Klebsiella pneumoniae*.
  - incubator
  - forceps
  - metric rulers
  - (70%) ethyl alcohol and beakers
  - Bunsen burner
- 

## **Method:**

### First Period

- 1) inoculate the surface of four Mueller-Hinton plates with *S. aureus*, *E. coli*, *P. aeruginosa*, and *K. pneumoniae*, respectively. Use a separate, sterile cotton swab for each bacterium. The swab is immersed in the culture tube, and the excess culture is squeezed on the inner side of the test tube.
- 2) The swab is then taken and streaked on the surface of the Mueller-Hinton plate three times, rotating the plate 60° after each streaking. Finally, run the swab around the edge of the agar.
- 3) Allow the culture to dry on the plate for 5 to 10 minutes at room temperature with the top in place.
- 3) Dispense the antibiotics onto the plate with alcohol-flamed Forceps.
- 4) gently pressing the disk with alcohol-flamed Forceps, make sure that contact is made between the antibiotic disk and the culture.
- 5) Incubate the plates for 16 to 18 hours at 35°C.

## Second Period

- 1) The diameter of the **inhibition zone** around the disk is measured to the nearest millimeter for each of the antibiotics tested, (by using metric rulers)
  - 2) Determine whether each organism is **susceptible, moderately susceptible, intermediate, or resistant to each chemotherapeutic agent** using the standardized table .
  - 3) record your results for every antibiotic on the table .
- 

Thank you ▶

